

Real-time monitoring of AAV Empty:Full Capsid Separations using Near infrared High-Precision Tunable Laser Spectroscopy (HPTLS)


Bryan Hassell, Ph.D., Marija Iloska, Ph.D., Jason Snyder, Hannah Furrelle, Larissa Miropolsky Nirrin Technologies



Motivation

Accurate measurement of AAV empty and full capsids is critical to gene therapy manufacturing, yet current tools-such as UV260/280, mass photometry, AUC, and TEM-are slow, offline, or insensitive to partially filled capsids1-3. UV absorbance methods, while fast, often fail to distinguish partial genomes from fully packaged vectors due to overlapping protein and nucleic acid signals².


We demonstrate the use of Near-Infrared High-Precision Tunable Laser Spectroscopy (NIR-HPTLS) to monitor AAV empty:full capsid separations in real time. By resolving subtle spectral differences between protein and nucleic acid content, our platform enables rapid, reagent-free fractionation analysis directly during chromatographic purification. This immediate feedback accelerates optimization, boosts yield, and opens the door to real-time release strategies.





> Nirrin makes the only combination band tunable laser capable of 200-300nm tuning range

High-Precision Tunable Laser Spectroscopy (HPTLS)



100X higher power density than broadband light source

> Yields wider dynamic range than other technologies Proprietary wavelength and amplitude referencing provides a calibrated spectrum with each scan

- > Allows for accurate library transfer across systems > This is the best-in-class signal-to-noise (SNR)
- spectrometer currently available (±10 µAU repeatability across spectrum)


Experimental setup

Äkta Go™ 11V 260/280 Nirrin HPTLS

1. Characterize "full" samples from AAV8 and AAV9 by injecting directly into Nirrin flow cell

Run AAV8/9 AEX with potassium chloride gradient

Monitor and compare in-line UV 260/280 (nucleic acid/protein) to HPTLS characterized signature



salt gradient, where increasing ionic strength displaces capsids based on surface charge. Traditional UV methods detect total absorbance but cannot resolve overlapping signals from salt, protein, and DNA. Using combination band near-infrared (NIR) spectroscopy, we simultaneously monitor the gradient profile (A, C) and extract precise spectral signatures for empty and full capsids (B) (Empty spectral information – E; Full spectral information - F).

Figure 2. Characterized "partial" run: Traditional analytical methods often misclassify partially filled AAV capsids

Pigure 2. octariatereeu partai fun: indiucina alayilay dina tenduo soten imisculos soten imisculos partai nue a Av Capata due to overlapito protein and DNA signals, espectandi yene using UU260/280 or bulk absorbance-based techniques. Partials may co-obiette with bland lagasids during RAX, complicating yield and potency assessments. NR HPTLS exploite Stintc combination and assorption features from nucleica add and potenci vibrational overtones (e.g., NH, CH, and OH stretches in the 4,200–5,000 cm⁻¹ region) to resolve subtle spectral differences between full and partially filled capsids124.5 (D) (Empty spectral information – E; Full spectral information – F; partial spectral information, P)

Conclusion

- NIR-HPTLS enables real-time monitoring of AAV empty:full separations
- Provides rapid, label-free analysis with no reagents or offline prep
- Detects subtle spectral differences between capsid protein and genome content.
- Enables fraction-level feedback during chromatography for improved control
- Supports process optimization and advances toward real-time release testing



1. Wright, J. F. (2020). Manufacturing and characterizing AAV-based vectors for use in clinical studies. Gene Therapy, 27(9), 317–327. Sommeregger, W. et al. (2022). A4V capsid content analysis using UV and orthogonal methods. Biotechnology and Bioengineering, 119(6), 1402–1413.
Metzsch, M. et al. (2020). Analytical ultracentrifugation and electron microscopy methods for A4V capsid characterization. Molecular Therapy – Method nent, 19, 447–457.

4. Wang, L. et al. (2018). Near-infrared spectroscopy of nucleic acids in gaueous solution. Vibrational Spectroscopy, 96, 1–9. 5. Huck, C. W. (2015). Near-infrared spectroscopy in bio-applications. Molecules, 20(3), 2587-2615.